

# Miltex

## HI-LIGHTS

### NO SCALPEL VASECTOMY INSTRUMENTS

#### Features

- The No Scalpel Vasectomy Ring Clamp isolates and secures the vas deferens without penetrating the skin
- The No Scalpel Vasectomy Hemostat pierces the scrotal sac to expose the vas deferens
- The vas deferens is lifted out of the scrotum with the No Scalpel Vasectomy Hemostat and occluded
- Gold finger rings for easy identification

*\*The No Scalpel Vasectomy technique is described in the article, "The No Scalpel Vasectomy", the Journal of Urology, February 1991, vol. 145, pages 341-344*



29-805  
No Scalpel Vasectomy Ring Clamp  
5 1/2", 3.5mm ID



29-850  
No Scalpel Vasectomy Hemostat  
5 1/2", Sharp point, smooth jaws, curved

**29-805 No Scalpel Vasectomy Ring Clamp**  
5 1/2", 3.5mm ID

**29-850 No Scalpel Vasectomy Hemostat**  
5 1/2", sharp point, smooth jaws, curved

**29-855 No Scalpel Vasectomy Instrument Set -**  
one of each above

#### **Miltex, Inc.**

Toll Free Customer Service: 866.854.8300  
To Fax Orders Toll Free: 866.854.8400  
www.miltex.com

