

Recommended Suture Sizes For Integra® Miltex® Tungsten Carbide Needle Holders

Reference	Description	Suture Sizes														
		Smallest												Largest		
		11/0	10/0	9/0	8/0	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2	3	4
8-3TC	Par NH Del Ser Jaw TC 5	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-	-
8-4TC	Derf NH TC 4-3/4	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-5TC	Par NH Del Ser Jaw TC TC 7	-	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-
8-6TC	Webster NH TC 5	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
8-7A-TC	Webster NH TC 4-3/4	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-
8-7TC	Webster NH TC 4-3/4	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
8-8A-TC	Halsey NH TC 5 Serr	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-8TC	Halsey NH TC 5 Smooth	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
8-14A-TC	Olsen-Hegar NH TC Smth	✓	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-
8-14TC	Olsen-Hegar NH TC Serr	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-15TC	Olsen-Hegar NH TC 5.5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-16TC	Olsen-Hegar NH TC 6.5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-17TC	Olsen-Hegar NH TC 7.5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-28TC	Mathieu NH TC 5.5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-30TC	Mathieu NH TC 6-1/4	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-40TC	Baumgartner NH TC 5.5	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-42TC	Mayo-Hegar NH TC 5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-44TC	Mayo-Hegar NH TC 6	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-46TC	Mayo-Hegar NH TC 7	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-48TC	Mayo-Hegar NH TC 8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-49TC	Mayo-Hegar NH TC 12	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-50A-TC	Crile-Wood NH TC 6 Smt	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
8-50TC	Crile-Wood NH TC 6 Ser	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-51TC	Crile-Wood NH TC 7	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-52TC	Crile-Wood NH TC 8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-53TC	Crile-Wood NH TC 9	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-55TC	Crile-Wood NH TC 10-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-57TC	Crile-Wood NH TC 12	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-62TC	Cooley NH X-Del TC 6	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-63TC	Cooley NH Xdel TC 6-3/4	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-65TC	Intracard NH TC 5 Del	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-66TC	Intracard NH TC 6 Del	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-67TC	Intracard NH TC 7 Del	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-69TC	Intracard NH TC 9 Del	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-81TC	New Orleans NH TC 7	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-82TC	Julian NH TC Str	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	-
8-83TC	Julian NH TC Cvd	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	-
8-96TC	Sarot NH TC 7-1/8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-100TC	Heaney NH TC 8-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-101TC	Ryder NH TC 5	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-102TC	Stratte NH TC 9	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-103TC	Ryder NH TC 6	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-104TC	Ryder NH TC 7	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-105TC	Ryder NH TC 8	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-106TC	Ryder NH TC 9	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-107TC	Ryder NH TC 10-1/4	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-114TC	Sarot NH TC 10-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-130TC	Masson NH TC 10-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-134TC	Wangensteen NH TC 10-3/4	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-706TC	Webster NH TC Quick-Release	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-
8-744TC	Mayo-Hegar NH TC Quick-Release	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-746TC	Mayo-Hegar NH TC Quick-Release	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-748TC	Mayo-Hegar NH TC Quick-Release	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-750TC	Crile-Wood NH TC Quick-Release	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-
8-751TC	Crile-Wood NH TC Quick-Release	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-752TC	Crile-Wood NH TC Quick-Release	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
8-753TC	Crile-Wood NH TC Quick-Release	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓
18-1816TC	Castroviejo NH TC Rnd	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
18-1828A-TC	Castroviejo NH TC Serr	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
18-1828TC	Castroviejo NH TC Smth	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
18-1832A-TC	Castroviejo NH TC Cvd Serr	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
18-1832TC	Castroviejo NH TC Cvd Smth	✓	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-
17-1090TC	Jacobson NH TC 7 Str	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	-
17-1091TC	Jacobson TC 7 Cvd	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	-
17-1092TC	Jacobson NH TC 8-1/2 Str	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	-

Recommended Suture Sizes For Integra® Miltex® Stainless Steel Needle Holders

Reference	Description	Suture Sizes															
		Smallest		9/0	8/0	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2	3	Largest	
8D-36	Gardner (Halsey) NH 5 Ser	-	-	-	-	-	-	-	-	✓	✓	✓	✓	-	-	-	
8-96	Sarot NH 7-1/8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-94	Crile NH 6	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-80	New Orleans NH 7	-	-	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	
8-8	Halsey NH 5	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-7	Webster NH 4-3/4 Del.	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-6	Webster NH 5	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-53	Crile-Wood NH 9	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-52	Crile-Wood NH 8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-51	Crile-Wood NH 7	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-50	Crile-Wood NH 6	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-48	Mayo-Hegar NH 8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-46	Mayo-Hegar NH 7	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-44	Mayo-Hegar NH 6	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-42	Mayo-Hegar NH 5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-40	Baumgartner NH 5.5	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-	
8-4	Derf NH 4-3/4	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-	
8-32	Mathieu NH 7.5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	-	-	
8-30	Mathieu NH 6-1/4	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	-	-	
8-28	Mathieu NH 5.5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	-	-	
8-26	Boynton NH 5 Del.	-	✓	✓	✓	-	-	-	-	-	-	-	-	-	-	-	
8-2	Collier NH 5	-	-	-	-	-	✓	✓	✓	✓	-	-	-	-	-	-	
8-17	Olsen-Hegar NH 7-1/4	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-16	Olsen-Hegar NH 6-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-15	Olsen-Hegar NH 5-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	-	-	
8-14A	Olsen-Hegar NH 4-3/4 Smth	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-	
8-14	Olsen-Hegar NH 4-3/4 Serr	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-	
8-136	Wangensteen NH 10-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-134	Wangensteen NH 10-3/4	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-130	Masson NH 10-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-122	Fino NH 10-1/2	-	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	
8-12	Neivert NH 5	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-114	Sarot NH 10-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-102	Stratte NH 9	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-100	Heaney NH 8-1/4	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
18-1833	Castroviejo NH Cvd W/Lock D	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	
18-1832	Castroviejo NH Cvd W/Lock	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	
18-1831	Castroviejo NH Str W/Lock D	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	
18-1828	Castroviejo NH Str W/Lock	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	
18-1824	Castroviejo NH Cvd Wo Lock	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	
18-1820	Castroviejo NH Str Wo Lock	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	-	
17-980	Micro NH 71/8 Str Jaw W/o	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	
17-1025	Micro NH 71/8 Cvd Jaw W/L	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	
17-1024	Micro NH 7-1/8 Str Jaw W/	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	
17-1020	Micro NH 51/8 Cvd Jaw W/L	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	
17-1010	Micro NH 71/8 Cvd Jaw Wo	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	
17-1000	Micro NH 51/8 Cvd Jaw Wo	-	-	-	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	

Recommended Suture Sizes For Integra® Miltex® Ceram-A-Grip® Needle Holders

Reference	Description	Suture Sizes															
		Smallest		9/0	8/0	7/0	6/0	5/0	4/0	3/0	2/0	0	1	2	3	Largest	
8-4CG	C-A-G Derf NH 4-3/4	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-	
8-6CG	C-A-G Webster NH 5	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-7CG	C-A-G Webster NH 4-3/4	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-8CG	C-A-G Halsey NH Smth	✓	✓	✓	✓	✓	✓	-	-	-	-	-	-	-	-	-	
8-15CG	C-A-G Olsen-Hegar 5-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	-	-	
8-16CG	C-A-G Olsen-Hegar 6-1/2	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	-	
8-40CG	C-A-G Baumgartner 5.5	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-	
8-42CG	C-A-G Mayo-Hegar NH 5	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-44CG	C-A-G Mayo-Hegar NH 6	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-46CG	C-A-G Mayo-Hegar NH 7	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-48CG	C-A-G Mayo-Hegar NH 8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	
8-50CG	C-A-G Crile-Wood NH 6	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-	
8-51CG	C-A-G Crile-Wood NH 7	-	-	-	-	-	✓	✓	✓	-	-	-	-	-	-	-	
8-52CG	C-A-G Crile-Wood NH 8	-	-	-	-	-	-	-	✓	✓	✓	✓	✓	✓	✓	✓	